




Bone tamps and impactors are designed to assist in positioning and securing bone grafts or implants during orthopedic procedures. Bone tamps provide controlled pressure to shape or compact graft material, while impactors deliver precise force to seat implants without causing damage. Manufactured from durable surgical-grade materials, these instruments offer reliability, balance, and ergonomic handling for accurate surgical execution.

# BONE TAMPS AND IMPACTORS

## BONE TAMPS

Cross-serrated tip, 6-1/4" (16 cm).



CODE #	DESCRIPTION
ML46893	2 mm diameter
ML46894	3 mm
ML46895	4 mm
ML46896	5 mm
ML46897	6 mm
ML46898	8 mm
ML46899	10 mm

## BONE IMPACTOR

Designed to tap bone graft or bone parts into place with minimal trauma to the bone. The compactor head is nylon with a grid surface to prevent slippage.




CODE #	DESCRIPTION
ML69980	Anodized aluminum handle with replaceable nylon caps, head diameter 20 mm, 6-1/2" (16.5 cm)

Autoclavable to 250°F (121°C)

## BONE CHIP PACKERS

Stainless steel handle, 9-7/8" (25 cm).




CODE #	DESCRIPTION
ML49750	3 mm x 16 mm
ML49751	2 mm x 9 mm
ML49752	6 mm x 9 mm
ML49753	6 mm x 19 mm
ML49754	1/4" (5 mm)

# IMPACTORS



## ENDOSCOPIC BONE TAMP

Straight, stainless steel handle, 17" (43 cm).



CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH
ML55676	3/8" (10 mm)	0.16"	12.59"	0.280" (7.1 mm)	4-1/4"

## LUMBAR AND CERVICAL IMPACTORS AND DRIVERS



CODE #	DESCRIPTION
ML55370	Impactor, standard offset, 6 mm x 12 mm, 8-1/2" (21.5 cm)
ML55371	Impactor, central canal, 7 mm x 16 mm, 8-1/2" (21.5 cm)
ML55372	Impactor, lateral narrow, 6 mm x 13 mm, 8-1/2" (21.5 cm)
ML55373	Impactor, micro offset, 4 mm x 9.5 mm, 7" (18 cm)
ML55374	Impactor, micro lateral, 4 mm x 10 mm, 7" (18 cm)
ML55375	Graft driver, lumbar, guarded, 10 mm x 13 mm, 8-1/2" (21.5 cm)
ML55376	Graft driver, anterior, guarded, 7 mm x 13 mm, 8-1/2" (21.5 cm)
ML55377	Graft driver, lumbar, unguarded, 10 mm x 13 mm, 8-1/2" (21.5 cm)
ML55378	Graft driver, anterior, unguarded, 7 mm x 13 mm, 8-1/2" (21.5 cm)