





# CARTILAGE AND TENDON INSTRUMENTS

Cartilage and tendon instruments are specially designed for delicate handling and precise manipulation of soft tissues during orthopedic procedures. Our instruments ensure minimal trauma while providing accurate control for cutting, grasping, or repositioning cartilage and tendons. Crafted from high-quality surgical steel, they support surgeons in performing meticulous repairs with confidence and care.

# COBB ELEVATORS

## COBB SPINAL ELEVATORS, HEX HANDLE



All stainless steel and hexagonal handle. Sharp, 9-1/2" (24 cm).



CODE #	DESCRIPTION
ML49340	1/2" (13 mm)
ML49342	3/4" (19 mm)
ML49344	1" (25.5 mm)
ML49346	1-1/4" (32 mm)

## COBB SPINAL ELEVATORS, ALUMINUM HANDLE

All stainless steel working end and anodized aluminum handle. Sharp, 11" (28 cm).





CODE #	DESCRIPTION
ML46413	3/8" (10 mm)
ML46414	1/2" (12 mm)
ML46416	3/4" (19 mm)
ML46418	1" (25 mm)
ML46420	1-1/4" (32 mm)

# ENDOSCOPIC ELEVATORS

## ENDOSCOPIC RIB ELEVATORS

Sharp, with phenolic handle




Phenolic handles autoclavable to 290°F (143°C)

CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML5564 0	Angled 45 degrees	9 mm radius	4"	8"	0.45" (11.5 mm)	17" (43 cm)
ML5564 4	Curved 90 degrees	9 mm radius	4"	8"	0.45" (11.5 mm)	17" (43 cm)

## ENDOSCOPIC COBB ELEVATORS

Phenolic handle.





Phenolic handles autoclavable to 290°F (143°C)

CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML55620	1/4" (6 mm)	2.7"	9.3"	0.285" (7.2 mm)	5"	17" (43 cm)
ML55610	1/4" (6 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)
ML55612	3/8" (10 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)
ML55614	0.49" (12.5 mm)	3-1/2"	10"	0.45" (11.5 mm)	6-1/2"	20" (51 cm)

# SPECIALTY ELEVATORS

## DARRACH ELEVATORS


Useful in shoulder surgery. Blunt tips with serrations can be used to elevate or retract soft tissue and muscle from the operative site.



CODE #	DESCRIPTION
ML46902	Small, 1/2" (12 mm) wide at tip, 10-1/4" (26 cm)
ML46903	Medium, 3/4" (19 mm) wide at tip, 10-1/4" (26 cm)
ML46904	Large, 1" (25 mm) wide at tip, 13-3/4" (35 cm)

## CHANDLER ELEVATORS


Blunt, slightly curved.



CODE #	TIP WIDTH	LENGTH
ML46250	8 mm	7-1/8" (18 cm)
ML46252	12 mm	8" (20.5 cm)
ML46254	17 mm	9" (23 cm)
ML46256	19 mm	10" (25.5 cm)



## HARRINGTON SPINAL ELEVATOR

All stainless steel



CODE #	DESCRIPTION
ML55798	Sharp, 3/4" (19 mm) wide tip, 12" (30.5 cm)

## SACHS NERVE ELEVATOR



CODE #	DESCRIPTION
ML57618	3 mm blunt, 7-1/2" (19 cm)