

CURETTES


Curettes are precision instruments used for scraping or debriding bone and soft tissue in orthopedic procedures. Designed with sharp, scoop-shaped tips, they allow for controlled removal of unwanted tissue while preserving surrounding structures. Manufactured from surgical-grade stainless steel, our curettes offer excellent balance, durability, and tactile feedback for precise surgical work.



BRUN CURETTES

BRUN (SPRATT) CURETTES

Oval cup, hollow handle




CODE #	TIP SHAPE	SIZE	LENGTH
ML46511	Straight	#5/0	7" (18 cm)
ML46512	Straight	#4/0	7" (18 cm)
ML46513	Straight	#3/0	7" (18 cm)
ML46515	Straight	#2/0	7" (18 cm)
ML46517	Straight	#0	7" (18 cm)
ML46519	Straight	#1	7" (18 cm)
ML46521	Straight	#2	7" (18 cm)
ML46523	Straight	#3	7" (18 cm)
ML46525	Straight	#4	7" (18 cm)
ML46529	Straight	#5	7" (18 cm)
ML46531	Straight	#6	7" (18 cm)
ML46563	Straight	#3/0	9" (23 cm)
ML46564	Straight	#2/0	9" (23 cm)
ML46565	Straight	#0	9" (23 cm)
ML46566	Straight	#1	9" (23 cm)
ML46567	Straight	#2	9" (23 cm)
ML46568	Straight	#3	9" (23 cm)
ML46569	Straight	#4	9" (23 cm)
ML46570	Straight	#5	9" (23 cm)
ML46571	Straight	#6	9" (23 cm)

CODE #	TIP SHAPE	SIZE	LENGTH
ML49110	Angled 45 degrees	#3/0	9" (23 cm)
ML49111	Angled 45 degrees	#2/0	9" (23 cm)
ML49112	Angled 45 degrees	#0	9" (23 cm)
ML49113	Angled 45 degrees	#1	9" (23 cm)
ML49114	Angled 45 degrees	#2	9" (23 cm)
ML49115	Angled 45 degrees	#3	9" (23 cm)
ML49116	Angled 45 degrees	#4	9" (23 cm)
ML49117	Angled 45 degrees	#5	9" (23 cm)
ML49118	Angled 45 degrees	#6	9" (23 cm)

BRUN (SPRATT) CURETTES





Round up, hollow handle



CODE #	TIP SHAPE	SIZE	LENGTH
ML49120	Straight	#3/0	6-3/4" (17 cm)
ML49121	Straight	#2/0	6-3/4" (17 cm)
ML49122	Straight	#0	6-3/4" (17 cm)
ML49123	Straight	#1	6-3/4" (17 cm)
ML49124	Straight	#2	6-3/4" (17 cm)
ML49125	Straight	#3	6-3/4" (17 cm)
ML49126	Straight	#4	6-3/4" (17 cm)

BRUN CURETTES


8-7/8" (22.5 cm).



CODE #	DESCRIPTION
343500	Size 000
343501	Size 00
343501	Size 0
343503	Size 1
343504	Size 2
343505	Size 3
343506	Size 4

BRUN CURETTES

Solid hexagonal handle, 9" (23 cm).



CODE #	TIP SHAPE	SIZE
ML46551	Straight	#3/0
ML46552	Straight	#2/0
ML46553	Straight	#0
ML46554	Straight	#1
ML46555	Straight	#2
ML46556	Straight	#3
ML46557	Straight	#4
ML46558	Straight	#5
ML46559	Straight	#6

CODE #	TIP SHAPE	SIZE
ML49150	Angled 45 degrees	#3/0
ML49151	Angled 45 degrees	#2/0
ML49152	Angled 45 degrees	#0
ML49153	Angled 45 degrees	#1
ML49154	Angled 45 degrees	#2
ML49155	Angled 45 degrees	#3
ML49156	Angled 45 degrees	#4
ML49157	Angled 45 degrees	#5
ML49158	Angled 45 degrees	#6

COBB CURETTES

COBB CURETTES, ALUMINUM HANDLE

All stainless steel working end and anodized aluminum, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE
ML46433	Angled	#3/0
ML46434	Angled	#2/0
ML46435	Angled	#0
ML46436	Angled	#1
ML46437	Angled	#2
ML46438	Angled	#3
ML46439	Angled	#4
ML46440	Angled	#5
ML46441	Angled	#6

CATALOG #	TIP SHAPE	SIZE
ML46422	Straight	#3/0
ML46423	Straight	#2/0
ML46424	Straight	#0
ML46425	Straight	#1
ML46426	Straight	#2
ML46427	Straight	#3
ML46428	Straight	#4
ML46429	Straight	#5
ML46430	Straight	#6

COBB CURETTES, STAINLESS STEEL HANDLE

All stainless steel working end and stainless steel, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE
ML46433S	Angled	#3/0
ML46434S	Angled	#2/0
ML46435S	Angled	#0
ML46436S	Angled	#1
ML46437S	Angled	#2
ML46438S	Angled	#3
ML46439S	Angled	#4
ML46440S	Angled	#5
ML46441S	Angled	#6

CODE #	TIP SHAPE	SIZE
ML46422S	Straight	#3/0
ML46423S	Straight	#2/0
ML46424S	Straight	#0
ML46425S	Straight	#1
ML46426S	Straight	#2
ML46427S	Straight	#3
ML46428S	Straight	#4
ML46429S	Straight	#5
ML46430S	Straight	#6

SPECIALTY BONE CURETTES


EPSTEIN CURETTES

Reverse angle, hollow handle, 8" (20.5 cm).




CODE #	DESCRIPTION
ML49580	#000
ML49582	#00
ML49584	#0

MARTINI BONE CURETTE



CODE #	DESCRIPTION
ML29620	Double-ended, 4 mm and 5 mm round cups, 5-1/2" (14 cm)

SCOVILLE CURETTES



CODE #	DESCRIPTION
ML49574	Straight, 10" (25.5 cm)
ML49576	Down angle, 10" (25.5 cm)
ML49578	Up angle, 10" (25.5 cm)

HATFIELD BONE CURETTES

Stainless steel with satin finish. Oval shaped tips, overall length 9" (23 cm). Complete set of nine curettes



CODE #	TIP SHAPE	SIZE	TIP LENGTH
584160	000	2.5 mm	4.0 mm
	Straight	2.8 mm	4.6 mm
	Straight	3.3 mm	5.0 mm
	Straight	3.8 mm	6.0 mm
	Straight	4.8 mm	7.6 mm
	Straight	5.8 mm	9.2 mm
	Straight	7.2 mm	11.6 mm
	Straight	8.7 mm	12.8 mm
	Straight	10.5 mm	15.5 mm


JANSEN CURETTE



CODE #	DESCRIPTION
ML29632	Double-ended, oval cups, 6" (15 cm)


VOLMLANN CURETTES

Double-ended, oval and round cups.



CODE #	DESCRIPTION
ML46560	5" (12.5 cm)
ML46561	8" (20.5 cm)
ML46562	6-1/2" (16.5 cm)


WILLIGER BONE CURETTE



CODE #	DESCRIPTION
ML29630	Double-ended, 3 mm and 4 mm oval cups, 5-1/2" (14 cm)

VOLMLANN CURETTES


Double-ended, oval cups.



CODE #	DESCRIPTION
ML49200	Double-ended, straight, 5" (12.5 cm)
ML49202	Double-ended, straight, 5 1/2" (14 cm)
ML46550	Double-ended, oval cups, 6 1/2" (16.5 cm)
ML49204	Double-ended, straight, 8" (20.5 cm)

RING CURETTES

Hollow handle with fenestrated cup.



CODE #	DESCRIPTION
ML55680	Size 0, small, 10" (25.5 cm)
ML55684	Size 0, small, 11-1/2" (29 cm)
ML55681	Size 2, medium, 10" (25.5 cm)
ML55685	Size 2, medium, 11-1/2" (29 cm)
ML55682	Size 6, large, 10-1/4" (26 cm)
ML55686	Size 6, large, 12" (30.5 cm)

CONE RING CURETTES

Knurled handle and fenestrated cup, 9" (23 cm).



CODE #	TIP SHAPE	SIZE	
ML49450	Straight	#1	8 mm diameter
ML49452	Straight	#2	6 mm diameter
ML49454	Straight	#3	3 mm diameter
ML49456	Angled 35 degrees	#1	8 mm diameter
ML49458	Angled 35 degrees	#2	6 mm diameter
ML49460	Angled 35 degrees	#3	3 mm diameter

RAY PITUITARY CURETTES

7-1/2" (19 cm).





CODE #	TIP SHAPE	SIZE
ML49403	Vertical angle 45 degrees	3 mm
ML49404	Vertical angle 45 degrees	4 mm
ML49405	Vertical angle 45 degrees	5 mm
ML49413	Vertical angle 90 degrees	3 mm
ML49414	Vertical angle 90 degrees	4 mm
ML49415	Vertical angle 90 degrees	5 mm
ML49423	Horizontal angle 45 degrees	3 mm
ML49424	Horizontal angle 45 degrees	4 mm
ML49425	Horizontal angle 45 degrees	5 mm
ML49433	Horizontal angle 90 degrees	3 mm
ML49434	Horizontal angle 90 degrees	4 mm
ML49435	Horizontal angle 90 degrees	5 mm

ENDOSCOPIC CURETTES


ENDO BRUN CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML55532	Angled 40°	2 x 4 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
ML55523	Straight	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
ML55533	Angled 40°	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
ML55543	Angled 40°	3 x 6 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)
ML55556	Angled 30°	6 x 9 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)


ENDO EPSTEIN CURETTES



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML55593	Reverseangled 40°	3 x 6 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
ML55586	Reverseangled 30°	6 x 9 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO SCOVILLE CURETTE


Endoscopic with phenolic handle.



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML55510	Straight	9 x 15 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO CONE RING CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
ML55573	Angled 40°	3 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
ML55566	Angled 30°	6 mm	10"	10"	0.45"	6-1/2"	20" (51 cm)

SMARTRACK CURETTE TRAYS

BRUN SMARTRACK FOR 7" CURETTES

Aluminum and stainless steel.



CODE #	DESCRIPTION
ML98160	Holds 12 7" MLeDedic® Brun Currettes

COBB SMARTRACK FOR 11" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
ML98170	Holds nine 11" MLeDedic® Cobb Currettes, Elevators or Gouges

HIBBS AND BRUN SMARTRACK FOR 9" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
ML98180	Holds 12 9" MLeDedic® Brun Currettes, or Hibbs Chisels, Gouges or Osteotomes